

# MATERIAL SPECIFICATIONS



MATERIAL	Manufacturer	Dk	Refractive Index	Wetting Angle	Specific Gravity	Shore D Hardness	Colors
SCLERAL LENS DESIGNS							
Tyro 97	Lagado	97	1.440	25	1.09	83	Clear, Light Blue
Optimum Extra	Contamac	100	1.431	3	1.16	75	Clear, Light Blue
Optimum Extreme	Contamac	125	1.432	6	1.15	77	Clear, Light Blue
CORNEAL CUSTOM GP DESIGNS							
SA-18	Lagado	18	1.490	25	1.12	87	Blue, Green
Optimum Classic	Contamac	26	1.450	12	1.19	83	Blue, Green, Gray
FLOSI	Lagado	26	1.430	23	1.27	86	Blue, Green, Gray, Brown
Hybrid FS	Contamac	31	1.446	N/A	1.17	79	Blue, Green, Gray
SA-32	Lagado	32	1.487	25	1.10	86	Blue, Green
ONSI-56	Lagado	56	1.452	7	1.20	86	Blue, Green, Gray, Brown
Optimum Comfort	Contamac	65	1.437	6	1.17	79	Blue, Green, Gray, Brown
Tyro-97	Lagado	97	1.440	25	1.09	83	Blue, Green, Gray, Brown
Optimum Extra	Contamac	100	1.431	3	1.16	75	Blue, Green, Clear, Gray
Optimum Extreme	Contamac	125	1.432	6	1.15	77	Blue, Green, Clear, Gray
CUSTOM							
PMMA	Lagado	0	1.495	25	1.19	93	Blue, Green, Gray, Brown